

AT A GLANCE

The EasyFoam® Kit consists of three parts:

- 10 ml low silicone disposable syringe filled with the required amount of sterile air with fixed bidirectional check valve and connector
- 5 ml low silicone disposable syringe (for injection purposes)
- needle 40 x 0.8 mm, 21G x 1½"



EasyFoam® Kit

Produce good foam easily



EasyFoam® Kit
for optimal
microfoam production

kreussler
PHARMA

EASYFOAM® KIT

ADVANTAGES

EasyFoam® is the syringe system for producing standardised microfoam. Using this system you can produce foam with the optimal ratio of gas to sclerosant in the mixture.

Foam sclerotherapy of varicose veins with Aethoxysklerol® (1% - 3%) is officially approved and recognised by the German Federal Institute for Drugs and Medical Devices (BfArM).

Double syringe system with fixed bidirectional check valve, connector and needle



- sterile packaging
- contains sterile gas for an optimal gas to liquid ratio in the mixture
- fixed bidirectional check valve
 - no leakage of sterile gas
- easy handling, low silicone syringes
 - silicone disrupts the foam
- optimal microfoam: standardised, stable, homogeneous, viscous and with fine bubbles
- long half-life of foam
- good usability proved in the marketing authorisation study (ESAF)

If you have any questions, please contact:

Chemische Fabrik Kreussler & Co. GmbH

P.O. Box 12 04 54 · 65082 Wiesbaden

Germany

Tel.: +49 (0)611 9271-0

Fax: +49 (0)611 9271-111

info@kreussler.com