


Aethoxysklerol®



Besenreiservarizen
Spider veins
Télangiectasies
Teleangectasia
Telangiectásias



Retikuläre und kleine Varizen
Reticular and small varices
Veines réticulaires et petites
Varici reticulares e piccolo
calibro
Varices reticulares y varices
pequeñas



Mittelgroße Varizen
Medium-sized varices
Varices moyennes
Varici di medio calibro
Varices medianas



Große Varizen
Large varices
Varices tronculaires
Varici di grosso calibro
Varices grandes




Hämorrhoiden
Haemorrhoids
Hémorroïdes
Emorroidi
Hemorroides

maximum number of ampoules to reach the maximum dose of 2 mg polidocanol per kg bodyweight in a patient weighing

50 kg 60 kg 75 kg 90 kg




2,5 mg/ml Polidocanol¹ 20 24 30 36



5 mg/ml Polidocanol¹ 10 12 15 18








10 mg/ml Polidocanol¹ 5 6 7,5 9

20 mg/ml Polidocanol¹ 2,5 3 3,5 4,5

30 mg/ml Polidocanol¹ 1,5 2 2,5 3

Volume (ml) per injection

0,1 – 0,2 ml

0,1 – 0,3 ml

0,5 – 2,0 ml³

0,5 – 2,0 ml³

– 1,0 ml⁴

Aethoxysklerol® universell einsetzbar bei optimaler Wirkung und Verträglichkeit
Aethoxysklerol® universally applicable with an optimum of efficacy and tolerance
Aetoxisclerol® universellement applicable avec efficacité et tolérance optimales
Atossisclerol® utilizzo universale con un massimo di efficacia e tolleranza
Etosisclerol® aplicación universal con efecto y tolerancia óptima

Availability of different concentrations depends on national marketing authorisation.

¹ Polidocanol, also known as Lauromacrogol 400, is produced by Kreussler and GMP-approved

² For spider veins, Aethoxysklerol 1% is administered only in central veins of spider veins

³ In the first treatment session, only one injection of 0,5 – 1 ml should be given

⁴ For sclerotherapy of haemorrhoidal disease, a total of 3 ml should not be exceeded during one treatment session; when treating an 11 o'clock haemorrhoid in men, the quantity injected must not exceed 0,5 ml